

# CNTD2 Antibody



PACO22035

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

2 isoforms of the human protein are produced by alternative splicing.

**Gene ID:**

CNTD2

**Uniprot**

Q9H8S5

**Synonyms:**

Cyclin N-terminal domain-containing protein 2; CNTD2;

**Immunogen:**

Synthesized peptide derived from internal of human CNTD2.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using CNTD2 antibody.

