

# COL23A1 Antibody



PACO22038

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:100-1:500

**Protein Background:**

COL23A1 is a member of the transmembrane collagens, a subfamily of the nonfibrillar collagens that contain a single pass hydrophobic transmembrane domain (Banyard et al., 2003). [supplied by OMIM, Mar 2008].

**Gene ID:**

COL23A1

**Uniprot**

Q86Y22

**Synonyms:**

Collagen alpha-1(XXIII) chain; COL23A1;

**Immunogen:**

Synthesized peptide derived from internal of human Collagen XXIII alpha 1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunofluorescence analysis of HepG2 cells, using Collagen XXIII alpha 1 antibody.

