

COX19 Antibody



PACO22039

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

May be required for the assembly of mitochondrial cytochrome c oxidase By similarity.

Gene ID:

COX19

Uniprot

Q49B96

Synonyms:

Cytochrome c oxidase assembly protein COX19; Hcox19; COX19;

Immunogen:

Synthesized peptide derived from internal of human COX19.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using COX19 antibody.

