

PACO22052

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Probable chaperone.

Gene ID:

DNAJB4

Uniprot

Q9UDY4

Synonyms:

DJB4; DnaJ homolog subfamily B member 4; DNAJW; DNJB4; Heat shock 40 kDa protein 1 homolog


Immunogen:

Synthesized peptide derived from internal of human DNAJB4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using DNAJB4 antibody.