

# RPS19BP1 Antibody



PACO22093

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

Direct regulator of SIRT1. Enhances SIRT1-mediated deacetylation of p53/TP53, thereby participating in inhibition of p53/TP53-mediated transcriptional activity.

**Gene ID:**

RPS19BP1

**Uniprot**

Q86WX3

**Synonyms:**

40S ribosomal protein S19-binding protein 1; RPS19-binding protein 1; S19BP;

**Immunogen:**

Synthesized peptide derived from internal of human RPS19BP1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human brain tissue using RPS19BP1 antibody.

