

ABHD14A Antibody



PACO22111

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Possible role in granule neuron development By similarity.

Gene ID:

ABHD14A

Uniprot

Q9BUJ0

Synonyms:

Abhydrolase domain-containing protein 14A; EC 3. -. -. -;


Immunogen:

Synthesized peptide derived from internal of human ABHD14A.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from RAW264.7 cells, using ABHD14A antibody.