

EPHX3 Antibody



PACO22114

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Belongs to the AB hydrolase superfamily.

Gene ID:

EPHX3

Uniprot

Q9H6B9

Synonyms:

abhydrolase domain containing 9; FLJ22408; MGC131519;


Immunogen:

Synthesized peptide derived from internal of human ABHD9.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from MCF-7 cells, using ABHD9 antibody.