

AKR1CL1 Antibody



PACO22134

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:500

Protein Background:

Belongs to the aldo/keto reductase family.

Gene ID:

AKR1C8P

Uniprot

Q5T2L2

Synonyms:

Aldo-keto reductase family 1 member C-like protein 1; EC 1.1.1. -;

Immunogen:

Synthesized peptide derived from internal of human AKR1CL1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunofluorescence analysis of A549 cells, using AKR1CL1 antibody.

