

AMPD2 Antibody



PACO22139

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

AMP deaminase plays a critical role in energy metabolism.

Gene ID:

AMPD2

Uniprot

Q01433

Synonyms:

AMD2; AMP deaminase 2; AMP deaminase isoform L;


Immunogen:

Synthesized peptide derived from internal of human AMPD2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, NIH-3T3 cells and COLO205 cells, using AMPD2 antibody.