

AIG1 Antibody



PACO22146

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play a role in androgen-regulated growth of hair follicles.

Gene ID:

AIG1

Uniprot

Q9NVV5

Synonyms:

Androgen-induced protein 1; AIG-1; CGI-103;


Immunogen:

Synthesized peptide derived from internal of human AIG1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and HeLa cells, using AIG1 antibody.