

ABCA6 Antibody



PACO22168

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Probable transporter which may play a role in macrophage lipid homeostasis.

Gene ID:

ABCA6

Uniprot

Q8N139

Synonyms:

ABC transporter ABCA6; ABCA6; ATP-binding cassette A6; ATP-binding cassette sub-family A member 6; ATP-binding cassette

Immunogen:

Synthesized peptide derived from C-terminal of human ABCA6.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, using ABCA6 antibody.