

DDX3Y Antibody



PACO22177

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Probable ATP-dependent RNA helicase. May play a role in spermatogenesis.

Gene ID:

DDX3Y

Uniprot

O15523

Synonyms:

ATP-dependent RNA helicase DDX3Y; DBY; DDXY; DEAD (Aspartate-glutamate-alanine-aspartate) box polypeptide; Y chromosome

Immunogen:

Synthesized peptide derived from internal of human DDX3Y.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using DDX3Y antibody.

