

BAGE2 Antibody



PACO22185

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Unknown. Candidate gene encoding tumor antigens.

Gene ID:

BAGE2

Uniprot

Q86Y30

Synonyms:

Cancer/testis antigen 2.2; CT2.2;

Immunogen:

Synthesized peptide derived from internal of human BAGE2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human testis tissue using BAGE2 antibody.

