

CA5A Antibody



PACO22234

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Reversible hydration of carbon dioxide. Low activity.

Gene ID:

CA5A

Uniprot

P35218

Synonyms:

ca5; ca5a; carbonic anhydrase va; mitochondrial; cav


Immunogen:

Synthesized peptide derived from internal of human CA5A.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A549 cells, LOVO cells and K562 cells, using CA5A antibody.