

CBLN3 Antibody



PACO22258

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

May be involved in synaptic functions in the CNS By similarity.

Gene ID:

CBLN3

Uniprot

Q6UW01

Synonyms:

Cerebellin-3;

Immunogen:

Synthesized peptide derived from internal of human CBLN3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using CBLN3 antibody.

