

C1QL2 Antibody



PACO22281

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May regulate the number of excitatory synapses that are formed on hippocampus neurons. Has no effect on inhibitory synapses By similarity.

Gene ID:

C1QL2

Uniprot

Q7Z5L3

Synonyms:

Complement C1q-like protein 2;


Immunogen:

Synthesized peptide derived from internal of human C1QL2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells and HeLa cells, using C1QL2 antibody.