

CPNE8 Antibody



PACO22288

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

May function in membrane trafficking. Exhibits calcium-dependent phospholipid binding properties By similarity.

Gene ID:

CPNE8

Uniprot

Q86YQ8

Synonyms:

copine VIII; copine-8;

Immunogen:

Synthesized peptide derived from N-terminal of human CPNE8.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using CPNE8 antibody.

