

PLA2G4E Antibody



PACO22299

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Calcium-dependent phospholipase A2 that selectively hydrolyzes glycerophospholipids in the sn-2 position By similarity.

Gene ID:

PLA2G4E

Uniprot

Q3MJ16

Synonyms:

Cytosolic phospholipase A2 epsilon; EC 3.1.1.4; cPLA2-epsilon; Phospholipase A2 group IVE;

Immunogen:

Synthesized peptide derived from internal of human PLA2G4E.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using PLA2G4E antibody.

