

# MARCH4 Antibody



PACO22319

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

**Gene ID:**

38047

**Uniprot**

Q9P2E8

**Synonyms:**

E3 ubiquitin-protein ligase MARCH4 ; EC 6.3.2. -; Membrane-associated RING finger protein 4; Membrane-associated RING-CH protein IV; MARCH-IV

**Immunogen:**

Synthesized peptide derived from internal of human MARCH4.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A