

PACO22324

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

Condensing enzyme that elongates saturated and monounsaturated very long chain fatty acids (VLCFAs). Highest activity toward C18 acyl-CoAs, especially C18:0 acyl-CoAs.

**Gene ID:**

ELOVL3

**Uniprot**

Q9HB03

**Synonyms:**

Elongation of very long chain fatty acids protein 3; Cold-inducible glycoprotein of 30 kDa;

**Immunogen:**

Synthesized peptide derived from internal of human ELOVL3.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human placenta tissue using ELOVL3 antibody.

