

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Belongs to the LDLR family. 3 isoforms of the human protein are produced by alternative splicing.

Gene ID:

LDLRAD1

Uniprot

Q5T700

Synonyms:

Low-density lipoprotein receptor class A domain-containing protein 1; LDLRAD1;


Immunogen:

Synthesized peptide derived from internal of human LDLRAD1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and HepG2 cells, using LDLRAD1 antibody.