

LDLRAD2 Antibody



PACO22393

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Belongs to the LDLR family.

Gene ID:

LDLRAD2

Uniprot

Q5SZI1

Synonyms:

ldlrad2; low density lipoprotein receptor class a domain containing 2; lrad2;


Immunogen:

Synthesized peptide derived from internal of human LDLRAD2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, LOVO cells and 3T3 using LDLRAD2 antibody.