

PACO22460

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play an important role in neural organelle transport, and in transduction of nerve signals or in nerve growth. May play a role in neurite extension. UniProtKB P97799.

Gene ID:

NRSN1

Uniprot

Q8IZ57

Synonyms:

neurensin 1; nrn1; p24; vesicular membrane protein p24; vmp


Immunogen:

Synthesized peptide derived from internal of human NRSN1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using NRSN1 antibody.