

# NPAS4 Antibody



PACO22467

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

Acts as a transcriptional activator in the presence of ARNT. Can activate the CME (CNS midline enhancer) element and the expression of the drebrin gene.

**Gene ID:**

NPAS4

**Uniprot**

Q8IUM7

**Synonyms:**

Neuronal PAS domain-containing protein 4; Neuronal PAS4; HLH-PAS transcription factor NXF;

**Immunogen:**

Synthesized peptide derived from internal of human NPAS4.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human brain tissue using NPAS4 antibody.

