

PACO22494

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane.

**Gene ID:**

PXMP2

**Uniprot**

Q9NR77

**Synonyms:**

Peroxisomal membrane protein 2; 22 kDa peroxisomal membrane protein;

**Immunogen:**

Synthesized peptide derived from internal of human PXMP2.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PXMP2 antibody.

