

RABL2A Antibody



PACO22539

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

RABL2A

Uniprot

Q9UBK7

Synonyms:

Rab-like protein 2A;


Immunogen:

Synthesized peptide derived from internal of human RABL2A.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HuvEc and COLO205 cells, using RABL2A antibody.