

PPP2R3B Antibody



PACO22601

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

Gene ID:

PPP2R3B

Uniprot

Q9Y5P8

Synonyms:

NY-REN-8; P2R3B; PP2A B; PP2A; subunit B


Immunogen:

Synthesized peptide derived from internal of human P2R3B.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells, using P2R3B antibody.