

EEF1G Antibody



PACO22605

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Probably plays a role in anchoring the complex to other cellular components.

Gene ID:

EEF1G

Uniprot

P26641

Synonyms:

eEF-1B gamma; EF-1-gamma; EEF1G; elongation factor 1-gamma; eukaryotic translation elongation factor 1 gamma


Immunogen:

Synthesized peptide derived from internal of human EEF1G.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using EF1G antibody.