

TCEAL3/TCEAL5/TCEAL6 Antibody



PACO22654

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

TCEAL3/TCEAL5/TCEAL6

Uniprot

Q969E4/Q5H9L2/Q6IPX3

Synonyms:

Transcription elongation factor A protein-like 3; TCEA-like protein 3; Transcription elongation factor S-II protein-like 3; TCEAL3; MSTP072

Immunogen:

Synthesized peptide derived from internal of human TCEAL3/5/6.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCEAL3/5/6 antibody.

