

USP42 Antibody



PACO22671

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Deubiquitinating enzyme which may play an important role during spermatogenesis By similarity.

Gene ID:

USP42

Uniprot

Q9H9J4

Synonyms:

deubiquitinating enzyme 42; EC 3.1.2.15; UB42; ubiquitin carboxyl-terminal hydrolase 42; ubiquitin thiolesterase 42


Immunogen:

Synthesized peptide derived from internal of human USP42.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HT-29 cells, using USP42 antibody.