

# ANKRD20A1 Antibody



PACO22712

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

**Gene ID:**

ANKRD20A1

**Uniprot**

Q5TYW2

**Synonyms:**

A20A1; ANKRD20A; ANKRD20A1; ankyrin repeat domain 20 family; member A1

**Immunogen:**


Synthesized peptide derived from internal of human A20A1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from 3T3 cells, using A20A1 antibody.