

ZNF691 Antibody



PACO22713

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF691

Uniprot

Q5VV52

Synonyms:

Zinc finger protein 691; ZNF691;

Immunogen:

Synthesized peptide derived from Internal of human ZNF691.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using ZNF691 antibody.