

ZNF575 Antibody



PACO22718

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF575

Uniprot

Q86XF7

Synonyms:

Zinc finger protein 575; ZNF575;


Immunogen:

Synthesized peptide derived from N-terminal of human ZNF575.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from LOVO cells and HUVEC cells, using ZNF575 antibody.