

# FGFRL1 Antibody



PACO22744

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Has a negative effect on cell proliferation By similarity.

**Gene ID:**

FGFRL1

**Uniprot**

Q8N441

**Synonyms:**

FGF homologous factor receptor; FGF receptor-like protein 1; FGFR-5; FGFR-like protein; FGFR5

**Immunogen:**

Synthesized peptide derived from Internal of human FGFRL1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A