

ZAR1 Antibody



PACO22765

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Essential for female fertility. May play a role in the oocyte-to-embryo transition By similarity.

Gene ID:

ZAR1

Uniprot

Q86SH2

Synonyms:

Oocyte-specific maternal effect factor; ZAR1; zygote arrest 1; Zygote arrest protein 1;


Immunogen:

Synthesized peptide derived from C-terminal of human ZAR1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, using ZAR1 antibody.