

TMEM101 Antibody



PACO22773

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May activate NF-kappa-B signaling pathways.

Gene ID:

TMEM101

Uniprot

Q96IK0

Synonyms:

Transmembrane protein 101; Putative NF-kappa-B-activating protein 130; TMEM101;


Immunogen:

Synthesized peptide derived from internal of human TMEM101.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HuvEc cells, using TMEM101 antibody.