

PACO22783

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Specifically catalyzes the 21-hydroxylation of steroids. Required for the adrenal synthesis of mineralocorticoids and glucocorticoids.

**Gene ID:**

CYP21A2

**Uniprot**

P08686

**Synonyms:**

21-hydroxylase; 21-OHase; CA21H; CAH1; CP21A

**Immunogen:**


Synthesized peptide derived from C-terminal of human CYP21A2.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from Jurkat/293 cells, using CYP21A2 antibody.