

# LY75 Antibody



PACO22818

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Acts as an endocytic receptor to direct captured antigens from the extracellular space to a specialized antigen-processing compartment By similarity. Causes reduced proliferation of B-lymphocytes. UniProtKB Q60767.

**Gene ID:**

LY75

**Uniprot**

O60449

**Synonyms:**

C-type lectin domain family 13 member B; CD205; CLEC13B; DEC-205; gp200-MR6

**Immunogen:**

Synthesized peptide derived from Internal of human LY75.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COLO205 cells, using LY75 antibody.