

# ORF4 Antibody



PACO22854

---

## Product Information

**Size:**

100ul

**Reactivity:**

Virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

plays a role in initiation of gene silencing at the early stage of viral infection.

**Gene ID:**

ORF4

**Uniprot**

P15184

**Synonyms:****Immunogen:**

Peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A