

# CG8889-RA Antibody



PACO22855

---

## Product Information

**Size:**

100ul

**Reactivity:**

Drosophila melanogaster

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

CG8889 gene product from transcript CG8889-RA.

**Gene ID:**

CG8889-RA

**Uniprot**

D3DMZ6

**Synonyms:****Immunogen:**

Peptide sequence around aa.178~182 (K-R-Q-T-L) derived from Fruit fly RE57611p.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A