

JAK3 (Ab-785) Antibody



PACO23042

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:**Gene ID:**

JAK3

Uniprot

P52333

Synonyms:

EC 2.7.10.2; JAK-3; JAK3; Janus kinase 3; L-JAK


Immunogen:

Synthesized non-phosphopeptide derived from human JAK3 around the phosphorylation site of tyrosine 785 (S-D-Y(p)-E-L).

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using JAK3 (Ab-785) antibody.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using JAK3 (Ab-785) antibody.

