

ACTL6A Antibody



PACO23186

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IF:1:100-1:500

Protein Background:

Gene ID:

ACTL6A

Uniprot

O96019

Synonyms:

ACTL6; ACTL6A; MGC5382; actin-like 6A;


Immunogen:

Synthesized peptide derived from internal of human ACL6A.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using ACL6A antibody.

Immunofluorescence analysis of HUVEC cells, using ACL6A antibody.

