

ADCK4 Antibody



PACO23276

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IF:1:100-1:500

Protein Background:

May play a role in CoQ10 (COQ10A and/or COQ10B) biosynthesis, which is required for podocyte migration.

Gene ID:

COQ8B

Uniprot

Q96D53

Synonyms:

EC 2.7.11. - ; ADCK4; Uncharacterized aarF domain-containing protein kinase 4;


Immunogen:

Synthesized peptide derived from internal of human ADCK4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using ADCK4 antibody.

Immunofluorescence analysis of A549 cells, using ADCK4 antibody.

