

PACO23295

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100,
IF:1:100-1:500

Protein Background:

May play a key role in the development of multiple organ systems and particularly in cardiac development By similarity.

Gene ID:

ICK

Uniprot

Q9UPZ9

Synonyms:

EC 2.7.11.22; ICK; Intestinal cell kinase; KIAA0936; LCK2

Immunogen:


Synthesized peptide derived from internal of human ICK.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ICK antibody.



Immunofluorescence analysis of HeLa cells, using ICK antibody.

