

NECAB3 Antibody



PACO23379

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IF:1:100-1:500

Protein Background:

Inhibits the interaction of APBA2 with beta-amyloid precursor protein (APP), and hence allows formation of beta-amyloid.

Gene ID:

NECAB3

Uniprot

Q96P71

Synonyms:

AB2BP; amyloid beta A4 protein-binding family A member 2-binding; APBA2BP;
dJ63M2.4; dJ63M2.5


Immunogen:

Synthesized peptide derived from internal of human NECAB3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO cells, HepG2 cells and HT-29 cells, using NECAB3 antibody.

Immunofluorescence analysis of A549 cells, using NECAB3 antibody.

