

ANKRD20A3 Antibody



PACO23381

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

Gene ID:

ANKRD20A3

Uniprot

Q5VUR7

Synonyms:

Ankyrin repeat domain-containing protein 20A3;


Immunogen:

Synthesized peptide derived from internal of human ANKRD20A3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 3T3 cells, using ANKRD20A3 antibody.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ANKRD20A3 antibody.

