

CCDC102B Antibody



PACO23419

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

2 isoforms of the human protein are produced by alternative splicing.

Gene ID:

CCDC102B

Uniprot

Q68D86

Synonyms:

Coiled-coil domain-containing protein 102B;


Immunogen:

Synthesized peptide derived from internal of human CCDC102B.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HuvEc cells and COLO cells, using CCDC102B antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using CCDC102B antibody.