

Phospho-TYK2 (Tyr1054) Antibody



PACO23860

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TYK2

Uniprot

P29597

Synonyms:

tyrosine kinase 2; JTK1;


Immunogen:

Peptide sequence around phosphorylation site of tyrosine 1054 (H-E-Y(p)-Y-R) derived from Human TYK2.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HT29 cells untreated (lane 1) or treated with Anisomycin (lane 2) using TYK2 (Phospho-Tyr1054) Antibody.