

Phospho-PPP1R1B (Thr75) Antibody



PACO23863

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:200

Protein Background:

Inhibitor of protein-phosphatase 1.

Gene ID:

PPP1R1B

Uniprot

Q9UD71

Synonyms:

IPPD; PPP1R1B;


Immunogen:

Peptide sequence around phosphorylation site of threonine 75 (A-Y-T(p)-P-P) derived from Human DARPP-32.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunofluorescence staining of methanol-fixed HeLa cells using DARPP-32(Phospho-Thr75) Antibody.