

Phospho-STK11 (Thr189) Antibody



PACO23949

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

STK11

Uniprot

Q15831

Synonyms:

PJS, LKB1, hLKB1


Immunogen:

Peptide sequence around phosphorylation site of Threonine 189(G-G-T(p)-L-K) derived from Human LKB1.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 3T3 cells treated with PMA using LKB1 (Phospho-Thr189) Antibody. The lane on the right is treated with the antigen-specific peptide.